

# Electrocardiogram (ECG) IMPC\_ECG\_002

- [Purpose](#)
- [Experimental Design](#)
- [Procedure](#)
- [Notes](#)
- [Parameters](#)
- [Metadata](#)

## Purpose

To provide a high throughput method to obtain Electrocardiograms in a conscious mouse.

## Experimental Design

- **Minimum number of animals** : 5-7M + 5-7F
- **Age at test**: Ideal age = 12 weeks  $\pm$ 3 days.
- **Sex**: We would expect the results of this test to show sexual dimorphism

## Procedure

1. The lead plates are to be snapped into place onto the top of the pre-amplifier tower. The covering is removed to reveal three gel coated pads surrounded by a sticking plate. The plate will need to be covered with the extra cover in the package.
2. Turn on the combined amplifier and the pre-amplifier tower.
3. Double click the icon ECG acquisition on the acquisition computer.
4. Open the ECG set up file (for default settings).
6. Place mouse on pad, lowering the Red Acrylic Cubby to surround the mouse on 3 sides discouraging escape.
7. Press Start.
8. After the desired acquisition time, (5-10 minutes) stop the reading. There will be one long reading.
9. Save the data.
10. For additional readings create a new session using the same settings as before.
11. When saving sections with good readings, highlight the selected area and then save.

## Notes

## Data Analysis

1. Open Emouse Analyses icon
2. Select ECG signals
3. Choose folder (all readings in folder will show)
4. Click PNN X (for mice: N-N > than 6 ms)
5. Choose file(s) by highlighting
6. Go
7. Bottom file is the corrected file

8. Red dots should be on peak of R waves, if image appears inverted click invert
9. Click Add, or minus if R waves are not marked with red dots or if too many are marked
  - L click to zoom in
  - R click to zoom out
10. 'What if?' button to remove unwanted sections
  - L click image (zooms in)
  - L click left boundary
  - L click right boundary
11. Options- click more if want to exclude more sections
12. Undo available
13. Go
14. Here can input animal data if desired
15. Save- For the first mouse in in group, hit save, a new results folder will be created within the folder with the mouse data. Then can click quick save or next.
16. For the rest of the mice in the series, can hit quick save at this point- saves in last selected file - will group all files together in same excel sheet.
17. Open Emouse Analyses icon
18. Select ECG signals
19. Choose folder (all readings in folder will show)
20. Click PNN X (for mice: N-N > than 6 ms)
21. Choose file(s) by highlighting
22. Go
23. Bottom file is the corrected file
24. Red dots should be on peak of R waves, if image appears inverted click invert
25. Click Add, or minus if R waves are not marked with red dots or if too many are marked
  - L click to zoom in
  - R click to zoom out
26. 'What if?' button to remove unwanted sections
  - L click image (zooms in)
  - L click left boundary
  - L click right boundary
27. Options- click more if want to exclude more sections
28. Undo available
29. Go
30. Here can input animal data if desired
31. Save- For the first mouse in in group, hit save, a new results folder will be created within the folder with the mouse data. Then can click quick save or next
32. For the rest of the mice in the series, can hit quick save at this point- saves in last selected file - will group all files together in same excel sheet

Examples of good readings







Figure C. Analysis phase, with the options to remove sections on the 'What if?' button below.



## Parameters

	Version	Type	Req. Upload	Req. Analysis	Annotation	Increment	Option	Ontology Options	Derived	Unit	Data Type
<a href="#">Number of signals</a> IMPC_ECG_001_001	1.2	simpleParameter	✓								FLOAT
<a href="#">HR</a> IMPC_ECG_002_001	1.1	simpleParameter	✓		✓					bpm	FLOAT
<a href="#">CV</a> IMPC_ECG_003_001	1.0	simpleParameter			✓					%	FLOAT
<a href="#">RR</a> IMPC_ECG_004_001	1.2	simpleParameter	✓		✓					ms	FLOAT
<a href="#">PQ</a> IMPC_ECG_005_001	1.0	simpleParameter			✓					ms	FLOAT
<a href="#">PR</a> IMPC_ECG_006_001	1.1	simpleParameter	✓		✓					ms	FLOAT
<a href="#">QRS</a> IMPC_ECG_007_001	1.2	simpleParameter	✓		✓					ms	FLOAT
<a href="#">ST</a> IMPC_ECG_008_001	1.0	simpleParameter	✓		✓					ms	FLOAT
<a href="#">QTc</a> IMPC_ECG_009_002	2.0	simpleParameter	✓		✓					ms	FLOAT
<a href="#">HRV</a> IMPC_ECG_010_001	1.0	simpleParameter			✓					bpm	FLOAT
<a href="#">QTc Dispersion</a> IMPC_ECG_011_001	1.0	simpleParameter			✓					ms	FLOAT
<a href="#">Mean SR amplitude</a> IMPC_ECG_012_001	1.1	simpleParameter								mV	FLOAT
<a href="#">Mean R amplitude</a> IMPC_ECG_013_001	1.1	simpleParameter								mV	FLOAT
<a href="#">rMSSD</a> IMPC_ECG_014_001	1.0	simpleParameter			✓					ms	FLOAT
<a href="#">pNN5(6&gt;ms)</a> IMPC_ECG_015_001	1.2	simpleParameter								%	FLOAT
<a href="#">Waveform Image</a> IMPC_ECG_025_001	1.0	seriesMediaParameter									IMAGE
<a href="#">Waveform Image Comment</a> IMPC_ECG_026_001	1.0	simpleParameter									TEXT

## Metadata

	Version	Type	Req. Upload	Req. Analysis	Annotation	Increment	Option	Ontology Options	Derived	Unit	Data Type
<a href="#">Equipment ID</a> IMPC_ECG_016_001	1.0	procedureMetadata	✓								TEXT
<a href="#">Equipment Manufacturer</a> IMPC_ECG_017_001	1.0	procedureMetadata	✓	✓			Mouse Specifics, Inc. AD Instruments World Precision Instruments				TEXT
<a href="#">Equipment Model</a> IMPC_ECG_018_001	1.0	procedureMetadata	✓	✓			ECGenie ML870/p ML826/FE132 Iso-DAM8A ECGenie + gel pads ML866				TEXT
<a href="#">Anesthetic</a> IMPC_ECG_019_001	1.0	procedureMetadata	✓	✓			Isoflurane Avertin Tribromoethanol No anesthesia				TEXT
<a href="#">Experimenter ID</a> IMPC_ECG_020_001	1.0	procedureMetadata	✓								TEXT
<a href="#">Noise level</a> IMPC_ECG_021_001	1.0	procedureMetadata									TEXT
<a href="#">Light level</a> IMPC_ECG_022_001	1.0	procedureMetadata									TEXT
<a href="#">Date equipment last calibrated</a> IMPC_ECG_023_001	1.1	procedureMetadata									DATE
<a href="#">Analysis Software</a> IMPC_ECG_024_001	1.0	procedureMetadata		✓			eMouse Matlab				TEXT