

ALTERNATIVE - Grip-Strength ALTESLIM_009_001

Parameters and Metadata

ALTERNATIVE - Forelimb grip strength measurement ALTESLIM_009_001_001 | v1.0
seriesParameter

Req. Analysis: false Req. Upload: true Is Annotated: false

Unit Measured: g

Description: ALTERNATIVE - Forelimb_grip_strength_measurement

Increments: 1, 2, 3,

ALTERNATIVE - Forelimb and hindlimb grip strength measurement ALTESLIM_009_001_002 | v1.0
seriesParameter

Req. Analysis: false Req. Upload: true Is Annotated: false

Unit Measured: g

Description: ALTERNATIVE - Forelimb_and_hindlimb_grip_strength_measurement

Increments: 1, 2, 3,

ALTERNATIVE - Body weight ALTESLIM_009_001_003 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: g

Description: ALTERNATIVE - Body_weight

ALTERNATIVE - General comments about the mouse ALTE

SLIM_009_001_004 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - General_comments_about_the_mouse

ALTERNATIVE - Forelimb grip strength measurement mean ALTESLIM_009_001_701 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: g

Description: ALTERNATIVE - Forelimb_grip_strength_measurement_mean

Derivation:

archived('sum_of_increments(ESLIM_009_001_001)/number_of_increments
(ESLIM_009_001_001)')

ALTERNATIVE - Forelimb and hindlimb grip strength measurement mean ALTESLIM_009_001_702 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: g

Description: ALTERNATIVE - Forelimb_and_hind_limb_grip_strength_measurement_mean

Derivation:

archived('sum_of_increments(ESLIM_009_001_002)/number_of_increments
(ESLIM_009_001_002)')

ALTERNATIVE - Forelimb grip strength normalised against body weight ALTESLIM_009_001_703 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - Forelimb_grip_strength_normalised_against_body_weight

Derivation:

archived('(sum_of_increments(ESLIM_009_001_001)/number_of_increments
(ESLIM_009_001_001))/ESLIM_009_001_003')

ALTERNATIVE - Forelimb and hindlimb grip strength normalised against body weight ALTESLIM_009_001_704 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description:

ALTERNATIVE - Forelimb_and_hind_limb_grip_strength_normalised_against_body_weight

Derivation:

archived('(sum_of_increments(ESLIM_009_001_002)/number_of_increments
(ESLIM_009_001_002))/ESLIM_009_001_003')

ALTERNATIVE - Equipment name ALTESLIM_009_001_801 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - Equipment_name

ALTERNATIVE - Equipment manufacturer ALTESLIM_009_001_802 | v1.0

procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - Equipment_manufacturer

ALTERNATIVE - Equipment model ALTESLIM_009_001_803 | v1.0

procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - Equipment_model

ALTERNATIVE - Grid model ALTESLIM_009_001_804 | v1.0

procedureMetadata

Req. Analysis: true

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - Grid_model
