

ALTERNATIVE - Rotarod ALTESLIM_010_001

Parameters and Metadata

ALTERNATIVE - Latency to fall ALTESLIM_010_001_001 | v1.0

seriesParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: s

Description: ALTERNATIVE - Latency_to_fall

Increments: 1, 2, 3,

ALTERNATIVE - Passive rotation ALTESLIM_010_001_002 | v1.0

seriesParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - Passive_rotation

Increments: 1, 2, 3,

Options: no, yes,

ALTERNATIVE - Body weight ALTESLIM_010_001_003 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: g

Description: ALTERNATIVE - Body_weight

ALTERNATIVE - General comments about the mouse ALTE
SLIM_010_001_004 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - General_comments_about_the_mouse

ALTERNATIVE - Latency to fall mean ALTESLIM_010_001_701 | v1.
0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: s

Description: ALTERNATIVE - Latency_to_fall_mean

Derivation:

archived('sum_of_increments(ESLIM_010_001_001)/number_of_increments
(ESLIM_010_001_001)')

ALTERNATIVE - Equipment name ALTESLIM_010_001_801 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - Equipment_name

ALTERNATIVE - Equipment manufacturer ALTESLIM_010_001_802 | v1.0

procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - Equipment_manufacturer

ALTERNATIVE - Equipment model ALTESLIM_010_001_803 | v1.0

procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - Equipment_model

ALTERNATIVE - Surface of the rod ALTESLIM_010_001_804 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - Surface_of_the_rod

ALTERNATIVE - Diameter of the rod ALTESLIM_010_001_805 | v1.0

procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Unit Measured: cm

Description: ALTERNATIVE - Diameter_of_the_rod

ALTERNATIVE - Acceleration mode ALTESLIM_010_001_806 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - Acceleration_mode

ALTERNATIVE - Number of mice on the rod per run ALTESLI M_010_001_807 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - Number_of_mice_on_the_rod_per_run

ALTERNATIVE - First inter-trial interval ALTESLIM_010_001_808 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: min

Description: ALTERNATIVE - First_inter_trial_interval

ALTERNATIVE - Second inter-trial interval ALTESLIM_010_001_809 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: min

Description: ALTERNATIVE - Second_inter_trial_interval
