

# ALTERNATIVE - Light-Dark Test ALTIMPC\_LDT\_001

## Parameters and Metadata

### ALTERNATIVE - Latency to first transition into dark ALTIMPC

C\_LDT\_006\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: s

Description: ALTERNATIVE - latency\_to\_first\_transition\_into\_dark

---

### ALTERNATIVE - Percent time in dark ALTIMPC\_LDT\_007\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: %

Description: ALTERNATIVE - percent\_time\_in\_dark

#### Derivation:

```
mul(div(sumOfIncrements('IMPC_LDT_004_001',1), add(sumOfIncrements('IMPC_LDT_004_001',1), sub(sumOfIncrements('IMPC_LDT_002_001',1), 'IMPC_LDT_006_001'))), 100)
```

---

## ALTERNATIVE - Percent time in light ALTIMPC\_LDT\_008\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: %

Description: ALTERNATIVE - percent\_time\_in\_light

### Derivation:

```
mul(div(sub(sumOfIncrements('IMPC_LDT_002_001',1), 'IMPC_LDT_006_001'), add
(sumOfIncrements('IMPC_LDT_004_001',1), sub(sumOfIncrements('IMPC_LDT_002_001',
1), 'IMPC_LDT_006_001'))), 100)
```

---

## ALTERNATIVE - Transition ALTIMPC\_LDT\_009\_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - transition

Options: No, Yes,

---

## ALTERNATIVE - Light level light chamber ALTIMPC\_LDT\_016\_001

| v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

**Unit Measured:** Lux

**Description:** ALTERNATIVE - light\_level\_light\_chamber

**Options:** 100 - 200, 150 - 200, 275 Lux,

---

## **ALTERNATIVE - Light level dark chamber** ALTIMPC\_LDT\_017\_001 | v1.0

[procedureMetadata](#)

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** false

**Unit Measured:** Lux

**Description:** ALTERNATIVE - light\_level\_dark\_chamber

**Options:** 2 - 5, 4 - 7,

---

## **ALTERNATIVE - Experimenter ID** ALTIMPC\_LDT\_018\_001 | v1.0

[procedureMetadata](#)

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** false

**Description:** ALTERNATIVE - experimenter\_id

---

## **ALTERNATIVE - Start time** ALTIMPC\_LDT\_019\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - start\_time

---

## **ALTERNATIVE - Arena ID** ALTIMPC\_LDT\_020\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - arena\_id

---

## **ALTERNATIVE - Arena version** ALTIMPC\_LDT\_021\_001 | v1.0

procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - arena\_version

Options: Custom made, MED-OFAS-RSU, SuperFlex ECC,

---

## **ALTERNATIVE - Equipment ID** ALTIMPC\_LDT\_022\_001 | v1.0

procedureMetadata

**Req. Analysis:** true

**Req. Upload:** true

**Is Annotated:** false

**Description:** ALTERNATIVE - equipment\_id

**Options:** Custom made, Dark Box Insert ENV-516, Light-Dark Insert, SuperFlex Dark Enclosure,

---

## **ALTERNATIVE - Equipment manufacturer** ALTIMPC\_LDT\_023\_001 | v1.0

[procedureMetadata](#)

**Req. Analysis:** true

**Req. Upload:** true

**Is Annotated:** false

**Description:** ALTERNATIVE - equipment\_manufacturer

**Options:** Accuscan, Custom made, Med Associates, Med Association Inc., Omnitech,

---

## **ALTERNATIVE - Equipment model** ALTIMPC\_LDT\_024\_001 | v1.0

[procedureMetadata](#)

**Req. Analysis:** true

**Req. Upload:** true

**Is Annotated:** false

**Description:** ALTERNATIVE - equipment\_model

**Options:** Activity Test Chamber ENV-515, Custom made, SuperFlex Environmental Control Chamber, Versamax,

---

## ALTERNATIVE - Disinfectant ALTIMPC\_LDT\_025\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - disinfectant

Options: 1% Virkon solution, 10% Nolvasan Solution, 50% alcohol, 50% Ethanol, 70% ethanol, 70% Isopropanol, Alcohol, Clidox, Distel 2%, Pursept-A,

---

## ALTERNATIVE - Starting position ALTIMPC\_LDT\_026\_001 | v1.0

procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - starting\_position

Options: Facing away from entry, Facing towards entry,

---

## ALTERNATIVE - Light-Dark box's walls opaque ALTIMPC\_LDT\_027\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - light\_dark\_box\_s\_walls\_opaque

Options: No, Yes,

---

## **ALTERNATIVE - Light-Dark box set-up** ALTIMPC\_LDT\_028\_001 | v1.0

procedureMetadata

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** false

**Description:** ALTERNATIVE - light\_dark\_box\_set\_up

**Options:** Cabinet enclosure, Standalone,

---

## **ALTERNATIVE - Light side arena dimensions (L x W x H, in mm)** ALTIMPC\_LDT\_029\_001 | v1.0

procedureMetadata

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** false

**Description:** ALTERNATIVE - light\_side\_arena\_dimensions\_l\_x\_w\_x\_h\_in\_mm

**Options:** 180 x 207 x 140, 250 x 190 x 250, 408 x 202 x 300, 431.8 x 212.725 x 203.2, 433 x 215 x 305,

---

## **ALTERNATIVE - Dark side arena dimensions (L x W x H, in mm)** ALTIMPC\_LDT\_030\_001 | v1.0

procedureMetadata

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** false

**Description:** ALTERNATIVE - dark\_side\_arena\_dimensions\_l\_x\_w\_x\_h\_in\_mm

**Options:** 150 x 190 x 250, 180 x 207 x 140, 396 x 192 x 98, 425 x 210 x 305,  
431.8 x 212.725 x 203.2,

---

## ALTERNATIVE - Entry way dimensions (L x W x H, in mm)

ALTIMPC\_LDT\_031\_001 | v1.0

procedureMetadata

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** false

**Description:** ALTERNATIVE - entry\_way\_dimensions\_l\_x\_w\_x\_h\_in\_mm

**Options:** 101.6 x 3.175 x 101.6, 102 x 5 x 32, 50 x 100 x 5, 90 x 3 x 115,

---

## ALTERNATIVE - Side changes ALTIMPC\_LDT\_001\_001 | v1.0

seriesParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** false

**Unit Measured:** count

**Description:** ALTERNATIVE - null

**Increments:** 5, 20, 15, 10,

---

## ALTERNATIVE - Light side time spent ALTIMPC\_LDT\_002\_001 | v1.

0

seriesParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: s

Description: ALTERNATIVE - null

Increments: 10, 15, 5, 20,

---

## ALTERNATIVE - Time mobile light side ALTIMPC\_LDT\_003\_001 | v1

.0

seriesParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: s

Description: ALTERNATIVE - null

Increments: 20, 5, 10, 15,

---

## ALTERNATIVE - Dark side time spent ALTIMPC\_LDT\_004\_001 | v1.

0

seriesParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

**Unit Measured:** s

**Description:** ALTERNATIVE - null

**Increments:** 5, 15, 20, 10,

---

**ALTERNATIVE - Time mobile dark side** ALTIMPC\_LDT\_005\_001 | v1  
.0

seriesParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** false

**Unit Measured:** s

**Description:** ALTERNATIVE - null

**Increments:** 5, 10, 15, 20,

---

**ALTERNATIVE - Light side distance travelled** ALTIMPC\_LDT\_01  
0\_001 | v1.0

seriesParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** false

**Unit Measured:** cm

**Description:** ALTERNATIVE - null

**Increments:** 10, 20, 5, 15,

---

## **ALTERNATIVE - Dark side distance travelled** ALTIMPC\_LDT\_01

1\_001 | v1.0

seriesParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: cm

Description: ALTERNATIVE - null

Increments: 5, 10, 15, 20,

---

## **ALTERNATIVE - Vertical activity (rearing)** ALTIMPC\_LDT\_012\_001

| v1.0

seriesParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: count

Description: ALTERNATIVE - null

Increments: 5, 20, 10, 15,

---

## **ALTERNATIVE - Horizontal activity** ALTIMPC\_LDT\_013\_001 | v1.0

seriesParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: count

Description: ALTERNATIVE - null

Increments: 20, 5, 10, 15,

---

## **ALTERNATIVE - Repetitive beam break ('stereotypy counts')** ALTIMPC\_LDT\_014\_001 | v1.0

seriesParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: count

Description: ALTERNATIVE - null

Increments: 10, 20, 15, 5,

---

## **ALTERNATIVE - Total distance travelled** ALTIMPC\_LDT\_015\_001 | v1.0

seriesParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: cm

Description: ALTERNATIVE - null

Increments: 10, 15, 20, 5,

---

## **ALTERNATIVE - Time spent in the light side in the first 5 minutes** ALTIMPC\_LDT\_032\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** false

**Description:** ALTERNATIVE - null

**Derivation:**

```
ifElse('IMPC_LDT_006_001'<300, sub(incrementValue('IMPC_LDT_002_001',  
1),'IMPC_LDT_006_001'), 300)
```

---

## **ALTERNATIVE - Test duration** ALTIMPC\_LDT\_033\_001 | v1.0

procedureMetadata

**Req. Analysis:** true

**Req. Upload:** true

**Is Annotated:** false

**Unit Measured:** min

**Description:** ALTERNATIVE -

**Options:** 10, 15, 20, 5,

---

## **ALTERNATIVE - Software version** ALTIMPC\_LDT\_034\_001 | v1.0

procedureMetadata

**Req. Analysis:** true

**Req. Upload:** true

**Is Annotated:** false

**Description:** ALTERNATIVE - null

**Options:** Activity Monitor V6.02, Activity Monitor V7.0.5.10, Fusion v6.5,  
Versamap Version 1.7,

---