

ALTERNATIVE - Rotarod ALTICSLA_ROT_001

Parameters and Metadata

ALTERNATIVE - Equipment manufacturer ALTICSLA_ROT_005_01 | v1.0

procedureMetadata

Req. Analysis: false Req. Upload: true Is Annotated: false

Description: ALTERNATIVE - null

Options: Bioseb,

ALTERNATIVE - Experimenter ID ALTICSLA_ROT_007_001 | v1.0

procedureMetadata

Req. Analysis: false Req. Upload: true Is Annotated: false

Description: ALTERNATIVE - null

ALTERNATIVE - Body weight ALTICSLA_ROT_004_001 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: false Is Annotated: false

Unit Measured: g

Description: ALTERNATIVE - null

ALTERNATIVE - Latency to fall ALTICSLA_ROT_001_001 | v1.1
seriesParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: s

Description: ALTERNATIVE - null

Increments: 1, 3, 2,

ALTERNATIVE - Acceleration mode ALTICSLA_ROT_010_001 | v1.0
procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - null

Options: starting speed 4-final 40,

ALTERNATIVE - Rod surface ALTICSLA_ROT_009_001 | v1.0
procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - null

Options: Foam rubber,

ALTERNATIVE - First intertrial interval ALTICSLA_ROT_011_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: min

Description: ALTERNATIVE - null

Options: 15,

ALTERNATIVE - Rod diameter ALTICSLA_ROT_008_001 | v1.1

procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Unit Measured: cm

Description: ALTERNATIVE - null

Options: 4.5,

ALTERNATIVE - Passive rotation ALTICSLA_ROT_003_001 | v1.3

seriesParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE -

Increments: 2, 3, 1,

Options: yes, no,

ALTERNATIVE - Latency to fall mean ALTICSLA_ROT_002_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: s

Description: ALTERNATIVE - null

Derivation: meanOfIncrements('ICSLA_ROT_001_001', 3)

ALTERNATIVE - Second intertrial interval ALTICSLA_ROT_012_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: min

Description: ALTERNATIVE - null

Options: 15,

ALTERNATIVE - Equipment model ALTICSLA_ROT_006_001 | v1.0
procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - null

Options: LE 8200,
