

ALTERNATIVE - Hot Plate ALTICSLA_HOT_001

Parameters and Metadata

ALTERNATIVE - Type of recording ALTICSLA_HOT_008_001 | v1.0

procedureMetadata

Req. Analysis: true Req. Upload: false Is Annotated: false

Description: ALTERNATIVE - null

Options: Manual,

ALTERNATIVE - Hot plate latency ALTICSLA_HOT_003_001 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: false

Unit Measured: s

Description: ALTERNATIVE -

ALTERNATIVE - Type of response ALTICSLA_HOT_004_001 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: false

Description: ALTERNATIVE -

Options: jump, lick, other, shake,

ALTERNATIVE - Equipment model ALTICSLA_HOT_006_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - null

Options: LE 7406,

ALTERNATIVE - Equipment manufacturer ALTICSLA_HOT_005_0

01 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - null

Options: Bioseb,

ALTERNATIVE - Experimenter ID ALTICSLA_HOT_001_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - null

ALTERNATIVE - Body weight ALTICSLA_HOT_002_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: g

Description: ALTERNATIVE - null

ALTERNATIVE - Plate temperature ALTICSLA_HOT_009_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: C

Description: ALTERNATIVE - null

Options: 52,
