

ALTERNATIVE - Rotarod ALTHMGULA_ROT_001

Parameters and Metadata

ALTERNATIVE - Latency to fall ALTHMGULA_ROT_001_001 | v1.0

seriesParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: s

Description: ALTERNATIVE - null

Increments: 1, 2, 3,

ALTERNATIVE - Latency to fall_Mean ALTHMGULA_ROT_002_001

| v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: s

Description: ALTERNATIVE - null

Derivation: meanOfIncrements('HMGULA_ROT_001_001', 3)

ALTERNATIVE - Passive rotation ALTHMGULA_ROT_003_001 | v1.0

seriesParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - null

Increments: 1, 2, 3,

Options: f, j, p,

ALTERNATIVE - Body weight ALTHMGULA_ROT_004_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: g

Description: ALTERNATIVE -

ALTERNATIVE - Comment ALTHMGULA_ROT_005_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - null

ALTERNATIVE - Animals per run ALTHMGULA_ROT_006_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - null

Options: 1, 2, 3, 4, 5,

ALTERNATIVE - Equipment ID ALTHMGULA_ROT_007_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - null

Options: 1019, 1069, 1070,

ALTERNATIVE - Equipment manufacturer ALTHMGULA_ROT_008_001 | v1.0

procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - null

Options: Bioseb,

ALTERNATIVE - Equipment model ALTHMGULA_ROT_009_001 | v1.

0

[procedureMetadata](#)

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE - null

Options: LE 8200,

ALTERNATIVE - Rod surface ALTHMGULA_ROT_010_001 | v1.0

[procedureMetadata](#)

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - null

Options: Rubber foam, Without rubber foam,

ALTERNATIVE - First intertrial interval ALTHMGULA_ROT_011_001

| v1.0

[procedureMetadata](#)

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: min

Description: ALTERNATIVE - null

Options: 15,

ALTERNATIVE - Second intertrial interval ALTHMGULA_ROT_012_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: min

Description: ALTERNATIVE - null

Options: 15,

ALTERNATIVE - Acceleration mode ALTHMGULA_ROT_013_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: rpm

Description: ALTERNATIVE - null

Options: 4-40,

ALTERNATIVE - Rod diameter ALTHMGULA_ROT_014_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: cm

Description: ALTERNATIVE - null

Options: 3, 5,
