

# ALTERNATIVE - Formalin Test ALTIMPC\_FOR\_001

## Parameters and Metadata

**ALTERNATIVE - Total duration of lick/bite behaviour** ALTIMPC\_FOR\_001\_001 | v1.0  
simpleParameter

Req. Analysis: false      Req. Upload: true      Is Annotated: false

Unit Measured: s

Description: ALTERNATIVE -

Derivation: sumOfIncrements('IMPC\_FOR\_020\_001', 1)

---

**ALTERNATIVE - Total number of lick/bite events** ALTIMPC\_FOR\_002\_001 | v1.0  
simpleParameter

Req. Analysis: false      Req. Upload: true      Is Annotated: false

Description: ALTERNATIVE -

---

**ALTERNATIVE - Total duration of drag/limp behaviour** ALTI

MPC\_FOR\_003\_001 | v1.0

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: false

Unit Measured: s

Description: ALTERNATIVE -

Derivation: sumOfIncrements('IMPC\_FOR\_021\_001', 1)

**ALTERNATIVE - Total number of drag/limp events** ALTIMPC\_

FOR\_004\_001 | v1.0

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: false

Description: ALTERNATIVE -

**ALTERNATIVE - Total duration of flinch behaviour** ALTIMPC\_

FOR\_005\_001 | v1.0

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: false

Unit Measured: s

Description: ALTERNATIVE -

**Derivation:** sumOfIncrements('IMPC\_FOR\_022\_001', 1)

---

**ALTERNATIVE - Total number of flinch events** ALTIMPC\_FOR\_006\_001 | v1.0

simpleParameter

**Req. Analysis:** false      **Req. Upload:** false      **Is Annotated:** false

**Description:** ALTERNATIVE -

---

**ALTERNATIVE - Anaesthetic** ALTIMPC\_FOR\_007\_001 | v1.0

procedureMetadata

**Req. Analysis:** false      **Req. Upload:** true      **Is Annotated:** false

**Description:** ALTERNATIVE -

**Options:** Sevoflurane, Isoflurane,

---

**ALTERNATIVE - Formalin Concentration** ALTIMPC\_FOR\_008\_001 | v1.0

procedureMetadata

**Req. Analysis:** true      **Req. Upload:** true      **Is Annotated:** false

**Description:** ALTERNATIVE -

**Options:** 5%, 2.5%,

---

**ALTERNATIVE - Amount of formalin injected** ALTIMPC\_FOR\_00

9\_001 | v1.0

procedureMetadata

**Req. Analysis:** true

**Req. Upload:** true

**Is Annotated:** false

**Unit Measured:** ul

**Description:** ALTERNATIVE -

**Options:** 20, 30,

---

**ALTERNATIVE - Syringe Gauge** ALTIMPC\_FOR\_010\_001 | v1.0

procedureMetadata

**Req. Analysis:** true

**Req. Upload:** true

**Is Annotated:** false

**Description:** ALTERNATIVE -

**Options:** 29, 30,

---

**ALTERNATIVE - Site of formalin injection** ALTIMPC\_FOR\_011\_001

| v1.0

procedureMetadata

**Req. Analysis:** true

**Req. Upload:** true

**Is Annotated:** false

**Description:** ALTERNATIVE -

**Options:** Plantar surface of left hind paw, Plantar surface of right hind paw,

---

## ALTERNATIVE - Arena Manufacturer ALTIMPC\_FOR\_012\_001 | v1.0

procedureMetadata

**Req. Analysis:** true

**Req. Upload:** true

**Is Annotated:** false

**Description:** ALTERNATIVE -

**Options:** In-house, IITC Life Science,

---

## ALTERNATIVE - Size of testing arena ALTIMPC\_FOR\_013\_001 | v1.0

procedureMetadata

**Req. Analysis:** true

**Req. Upload:** true

**Is Annotated:** false

**Description:** ALTERNATIVE -

**Options:** 13cm (H) x 36cm (W) x 40 cm(L), 12.7 cm (H) x 21.9 cm (W) x 21.59 cm (L), 36cm (L) X 40cm (W) X 13cm (H),

---

**ALTERNATIVE - Recording Software** ALTIMPC\_FOR\_014\_001 | v1.0

procedureMetadata

**Req. Analysis:** true

**Req. Upload:** true

**Is Annotated:** false

**Description:** ALTERNATIVE -

**Options:** Home Cage Analyzer, Actual Analytics, Noldus Media Recorder v4,

---

**ALTERNATIVE - Data scoring method** ALTIMPC\_FOR\_015\_001 | v1

.0

procedureMetadata

**Req. Analysis:** true

**Req. Upload:** true

**Is Annotated:** false

**Description:** ALTERNATIVE -

**Options:** Human scorer, from video, Neural network, from video,

---

**ALTERNATIVE - Annotation Software** ALTIMPC\_FOR\_016\_001 | v1.

0

procedureMetadata

**Req. Analysis:** true

**Req. Upload:** true

**Is Annotated:** false

**Description:** ALTERNATIVE -

**Options:** Simple Video Coder ,  
Behavioral Observation Research Interactive Software (BORIS), In-house developed system,

---

## **ALTERNATIVE - Duration of observation period** ALTIMPC\_FO

R\_017\_001 | v1.0

procedureMetadata

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** false

**Unit Measured:** min

**Description:** ALTERNATIVE -

**Options:** 60 minutes, 90 minutes,

---

## **ALTERNATIVE - Camera Position** ALTIMPC\_FOR\_018\_001 | v1.0

procedureMetadata

**Req. Analysis:** true

**Req. Upload:** true

**Is Annotated:** false

**Description:** ALTERNATIVE -

**Options:** On the side, Below the arena,

---

## **ALTERNATIVE - Disinfectant** ALTIMPC\_FOR\_019\_001 | v1.0

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE -

Options: 35% Isopropanol, Virkon 1% & 70% Ethanol, 2% Distel,

---

## ALTERNATIVE - Duration of lick/bite behaviour ALTIMPC\_FOR \_020\_001 | v1.0

seriesParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: s

Description: ALTERNATIVE -

Increments: Minimum 1

---

## ALTERNATIVE - Duration of drag/limp behaviour ALTIMPC\_F OR\_021\_001 | v1.0

seriesParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: s

Description: ALTERNATIVE -



Increments: Minimum 1

---

**ALTERNATIVE - Duration of flinch behaviour** ALTIMPC\_FOR\_022\_001 | v1.0

seriesParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: false

Unit Measured: s

Description: ALTERNATIVE -

Increments: Minimum 1

---

**ALTERNATIVE - Duration of licking/biting behavior (10-60 minutes)** ALTIMPC\_FOR\_023\_001 | v1.0

simpleParameter

Req. Analysis: false      Req. Upload: false      Is Annotated: false

Unit Measured: s

Description: ALTERNATIVE -

**Derivation:**  
sum(incrementValue('IMPC\_FOR\_020\_001', '10'),incrementValue('IMPC\_FOR\_020\_001', '15'),incrementValue('IMPC\_FOR\_020\_001', '20'),incrementValue('IMPC\_FOR\_020\_001', '25'),incrementValue('IMPC\_FOR\_020\_001', '30'),incrementValue('IMPC\_FOR\_020\_001', '35'),incrementValue('IMPC\_FOR\_020\_001', '40'),incrementValue('IMPC\_FOR\_020\_001', '45'),incrementValue('IMPC\_FOR\_020\_001', '50'),incrementValue('IMPC\_FOR\_020\_001', '55'))

---

## ALTERNATIVE - Time per bin/increment ALTIMPC\_FOR\_024\_001 | v1.0

procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE -

Options: 5 minutes,

---

## ALTERNATIVE - Time bins excluded ALTIMPC\_FOR\_025\_001 | v1.0

procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Description: ALTERNATIVE -

Options: Bin/increment 0 (anaesthesia effect), None,

---

## ALTERNATIVE - Experimenter ID ALTIMPC\_FOR\_026\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE -

-----