

# MicroCT E9.5 IMPC\_EML\_002

## Purpose

To assess morphological defects in E9.5 embryos from lethal strains with MicroCT.

## Experimental Design

- Set up timed matings with heterozygous mice
- Day 0 is defined as the midpoint of the prior dark cycle following the identification of a copulation plug.
- Collect embryos at E9.5 with 3 minimum homozygous and 1 littermate control
- Ensure embryos are bleed as much as possible.
- Reconstruct the scans according to manufacturer's instructions.

## Procedure

1. **Set up timed mating with heterozygous animals. Dissect and collect  $\geq 3$  homozygote embryos.**
2. **Process the embryos using the centre's staining and embedding protocol . Centre-specific methods are recorded within a document**
3. **Scan and generate reconstructions of embryos.**

## Notes

Embryos may be processed for Histopathology.

## Parameters and Metadata

### Camera to Source IMPC\_EML\_010\_001 | v1.0

[procedureMetadata](#)

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: mm

---

### Equipment Version IMPC\_EML\_006\_001 | v1.0

[procedureMetadata](#)

Req. Analysis: false

Req. Upload: true

Is Annotated: false

---

## Image Pixel Size IMPC\_EML\_007\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured:  $\mu\text{m}$

---

## Date of Scan IMPC\_EML\_014\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

---

## Equipment Manufacturer IMPC\_EML\_004\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

---

## Equipment ID IMPC\_EML\_003\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

---

## E9.5 MicroCT Preparation protocol IMPC\_EML\_008\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

---

## Screenshots IMPC\_EML\_017\_001 | v1.1

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## Object to Source IMPC\_EML\_009\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: mm

---

## Source Current IMPC\_EML\_016\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: uA

---

## Date of Sacrifice

IMPC\_EML\_012\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

---

## Date Preparation Finished

IMPC\_EML\_013\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

---

## Embryo reconstruction

IMPC\_EML\_001\_001 | v1.2

seriesMediaParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

---

## Experimenter ID

IMPC\_EML\_002\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

---

## Source Voltage

IMPC\_EML\_015\_001 | v1.1

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: kV

---

## Date equipment last calibrated IMPC\_EML\_011\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## Equipment Model IMPC\_EML\_005\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

---

## Hindbrain IMPC\_EML\_018\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Options: present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## 4th ventricle IMPC\_EML\_019\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## **Mesencephalic vesicle** IMPC\_EML\_020\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## **Forebrain** IMPC\_EML\_021\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## **Telecenphalic (lateral) ventricle** IMPC\_EML\_022\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## Rhombomeres IMPC\_EML\_023\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Options: present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## Notochord IMPC\_EML\_024\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## Infundibular recess of 3rd ventricle IMPC\_EML\_025\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## Rathke's pouch IMPC\_EML\_026\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

**Options:** present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## **Optic vesicle** IMPC\_EML\_027\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## **Otocyst** IMPC\_EML\_028\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## **1st branchial arch** IMPC\_EML\_029\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** present, absent, abnormal, unobservable, ambiguous, imageOnly,

---



## 1st arch artery IMPC\_EML\_030\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## 2nd branchial arch IMPC\_EML\_031\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Options: present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## 2nd arch artery IMPC\_EML\_032\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## Oropharynx IMPC\_EML\_033\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Options: present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## 3rd branchial arch IMPC\_EML\_034\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## 3rd arch artery IMPC\_EML\_035\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Options:** present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## Aortic sac IMPC\_EML\_036\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## Common atrial chamber IMPC\_EML\_037\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## **Atrioventricular canal** IMPC\_EML\_038\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## **Heart ventricle** IMPC\_EML\_039\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## **Primitive heart tube cardiac muscle** IMPC\_EML\_040\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Options:** present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## Dorsal aorta IMPC\_EML\_041\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Options: present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## Hepatic primordium IMPC\_EML\_042\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## Hindgut IMPC\_EML\_043\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Options: present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## Forelimb bud IMPC\_EML\_044\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

**Options:** present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## **Neural tube** IMPC\_EML\_045\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## **Caudal neuropore** IMPC\_EML\_046\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## **Somite** IMPC\_EML\_047\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Options:** present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## Somite count IMPC\_EML\_048\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## Delayed embryonic development IMPC\_EML\_049\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Options: yes, no,

---

## Vitelline artery IMPC\_EML\_050\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## Umbilical vein IMPC\_EML\_051\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## Foregut IMPC\_EML\_052\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## Midgut IMPC\_EML\_053\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Options: present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## Pericardial cavity IMPC\_EML\_054\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Options: present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## Peritoneal cavity IMPC\_EML\_055\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

**Options:** present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## **Sinus venosus** IMPC\_EML\_056\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** true

**Options:** present, absent, abnormal, unobservable, ambiguous, imageOnly,

---

## **Annotator ID** IMPC\_EML\_057\_001 | v1.0

procedureMetadata

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** false

## **Date of Annotation** IMPC\_EML\_058\_001 | v1.0

procedureMetadata

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** false

---