Body Weight HMGULA_BWT_001

Purpose

The body weight test measures the weight of the mouse in a time series, allowing monitoring of its evolution; also, it is required in many other procedures.

Equipment

Laboratory Balance

Procedure

- Place a cage containing mice on an operation table
- Remove a mouse from its cage. Determine and record its weight.
- Return the mouse to the cage or continue with a procedure.

Notes

The body weight of each mouse is to be measured monthly. Optional additional weighs may also be recorded.

Parameters and Metadata

Body weight HMGULA_BWT_001_001 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: false

Unit Measured: g

General comments about the mouse HMGULA_BWT_002_001 | v1.0

Req. Analysis: false	Req. Upload: false	Is Annotated: false	
Equipment ID HMGU procedureMetadata	LA_BWT_003_001 v1.0		
Req. Analysis: false	Req. Upload: true	Is Annotated: false	
Equipment manufacturer HMGULA_BWT_004_001 v1.0 procedureMetadata			
Req. Analysis: false	Req. Upload: true	Is Annotated: false	
Options: Ohaus, Denver Instrument, Scientech, Mettler Toledo, Kern & Sohn GmbH, Sartorius, Acculab, Radwag,			
Experimenter ID HMGULA_BWT_005_001 v1.0 procedureMetadata			
Req. Analysis: false	Req. Upload: true	Is Annotated: false	

Req. Analysis: false	Req. Upload: false	Is Annotated: false		
Equipment model HMGULA_BWT_007_001 v1.0 procedureMetadata				
Req. Analysis: false	Req. Upload: true	Is Annotated: false		
Options: AV2101, SI-2002, SL-3100D, TP-202, PG3001-S, Adventurer Pro, 440-47N,				

Options: AV2101, SI-2002, SL-3100D, TP-202, PG3001-S, Adventurer Pro, 440-47N, AM100, XS802S, AB104-S, AV212C, PCB2000-1, 440-33, 440-47, 572-35, 770-14, EMB220-1, EMB500-1, EMB600-2, EW600-2M, BL310, Navigator 34120, AV213C, VIC-511, VIC-123, BP6100, CPA3202S, EMB 200-2, Entris 2201-1S, 201-10, Cubis,

.....

Body weight curve HMGULA_BWT_008_001 | v1.1

seriesMediaParameter

Req. Analysis: false Req. Upload: false Is Annotated: false

Unit Measured: g

Derivation: unimplemented(")

Increments: Minimum 1