

ALTERNATIVE - Spontaneous breathing (GMC) ALTGMC_920_001

Parameters and Metadata

ALTERNATIVE - Body weight ALTGMC_920_001_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: g

Description: ALTERNATIVE - Body_weight

ALTERNATIVE - age ALTGMC_920_001_002 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: days

Description: ALTERNATIVE - age

ALTERNATIVE - f-active ALTGMC_920_001_003 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: 1/min

Description: ALTERNATIVE - factive

ALTERNATIVE - TV-active ALTGMC_920_001_004 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: ml

Description: ALTERNATIVE - TVactive

ALTERNATIVE - MV-active ALTGMC_920_001_005 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: ml/min

Description: ALTERNATIVE - MVactive

ALTERNATIVE - PIF-active ALTGMC_920_001_006 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: ml/s

Description: ALTERNATIVE - PIFactive

ALTERNATIVE - PEF-active ALTGMC_920_001_007 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: ml/s

Description: ALTERNATIVE - PEFactive

ALTERNATIVE - Ti-active ALTGMC_920_001_008 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: ms

Description: ALTERNATIVE - Tiactive

ALTERNATIVE - Te-active ALTGMC_920_001_009 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: ms

Description: ALTERNATIVE - Teactive

ALTERNATIVE - Ti/TT-active ALTGMC_920_001_010 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - TiTTactive

ALTERNATIVE - MIF-active ALTGMC_920_001_011 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: ml/s

Description: ALTERNATIVE - MIFactive

ALTERNATIVE - MEF-active ALTGMC_920_001_012 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: ml/s

Description: ALTERNATIVE - MEFactive

ALTERNATIVE - sTV-active ALTGMC_920_001_013 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: ml/g

Description: ALTERNATIVE - sTVactive

ALTERNATIVE - sMV-active ALTGMC_920_001_014 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: ml/(min*g)

Description: ALTERNATIVE - sMVactive

ALTERNATIVE - mean-f ALTGMC_920_001_015 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: 1/min

Description: ALTERNATIVE - meanf

ALTERNATIVE - f-rest ALTGMC_920_001_016 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: 1/min

Description: ALTERNATIVE - frest

ALTERNATIVE - TV-rest ALTGMC_920_001_017 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: ml

Description: ALTERNATIVE - TVrest

ALTERNATIVE - MV-rest ALTGMC_920_001_018 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: ml/min

Description: ALTERNATIVE - MVrest

ALTERNATIVE - PIF-rest ALTGMC_920_001_019 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: ml/s

Description: ALTERNATIVE - PIFrest

ALTERNATIVE - PEF-rest ALTGMC_920_001_020 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: ml/s

Description: ALTERNATIVE - PEFrest

ALTERNATIVE - Ti-rest ALTGMC_920_001_021 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: false Is Annotated: false

Unit Measured: ms

Description: ALTERNATIVE - Tirest

ALTERNATIVE - Te-rest ALTGMC_920_001_022 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: false Is Annotated: false

Unit Measured: ms

Description: ALTERNATIVE - Terest

ALTERNATIVE - Ti/TT-rest ALTGMC_920_001_023 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: false Is Annotated: false

Description: ALTERNATIVE - TiTTrest

ALTERNATIVE - MIF-rest ALTGMC_920_001_024 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: ml/s

Description: ALTERNATIVE - MIFrest

ALTERNATIVE - MEF-rest ALTGMC_920_001_025 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: ml/s

Description: ALTERNATIVE - MEFrest

ALTERNATIVE - sTV-rest ALTGMC_920_001_026 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: ml/g

Description: ALTERNATIVE - sTVrest

ALTERNATIVE - sMV-rest ALTGMC_920_001_027 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: ml/(min*g)

Description: ALTERNATIVE - sMVrest

ALTERNATIVE - f-sleep ALTGMC_920_001_028 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: 1/min

Description: ALTERNATIVE - fsleep

ALTERNATIVE - TV-sleep ALTGMC_920_001_029 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: ml

Description: ALTERNATIVE - TVsleep

ALTERNATIVE - MV-sleep ALTGMC_920_001_030 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: ml/min

Description: ALTERNATIVE - MVsleep

ALTERNATIVE - PIF-sleep ALTGMC_920_001_031 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: ml/s

Description: ALTERNATIVE - PIFsleep

ALTERNATIVE - PEF-sleep ALTGMC_920_001_032 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: ml/s

Description: ALTERNATIVE - PEFsleep

ALTERNATIVE - Ti-sleep ALTGMC_920_001_033 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: ms

Description: ALTERNATIVE - Tisleep

ALTERNATIVE - Te-sleep ALTGMC_920_001_034 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: ms

Description: ALTERNATIVE - Tesleep

ALTERNATIVE - Ti/TT-sleep ALTGMC_920_001_035 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Description: ALTERNATIVE - TiTTsleep

ALTERNATIVE - MIF-sleep ALTGMC_920_001_036 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: ml/s

Description: ALTERNATIVE - MIFsleep

ALTERNATIVE - MEF-sleep ALTGMC_920_001_037 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: ml/s

Description: ALTERNATIVE - MEFsleep

ALTERNATIVE - sTV-sleep ALTGMC_920_001_038 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: ml/g

Description: ALTERNATIVE - sTVsleep

ALTERNATIVE - sMV-sleep ALTGMC_920_001_039 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** false **Is Annotated:** false

Unit Measured: ml/(min*g)

Description: ALTERNATIVE - sMVsleep
