ALTERNATIVE - Modified SHIRPA ALTHAROLD_042_001

Parameters and Metadata

ALTERNATIVE - startle response to click ALTHAROLD_042_001_001 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - startle_response_to click

Options: 0, 1, 2,

ALTERNATIVE - Tremor ALTHAROLD_042_001_002 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - Tremor

Options: 0, 1, 2, 3,

ALTERNATIVE - Defecation ALTHAROLD_042_001_003 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - Defecation

Options: 0, 1,

ALTERNATIVE - Urination ALTHAROLD_042_001_004 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - Urination

Options: 0, 1,

ALTERNATIVE - Palpebral closure ALTHAROLD_042_001_005 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - Palpebral_Closure

Options: 0, 1, 2,

ALTERNATIVE - Lacrimation ALTHAROLD_042_001_006 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** false **Is Annotated:** false

Description: ALTERNATIVE - Lacrimation

Options: 0, 1,

ALTERNATIVE - Transfer arousal ALTHAROLD_042_001_007 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - Transfer_Arousal

Options: 0, 1, 2, 3, 4, 5, 6,

ALTERNATIVE - Locomotor activity ALTHAROLD_042_001_008 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - Locomotor_activity

ALTERNATIVE - Gait ALTHAROLD_042_001_009 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - Gait

Options: 0, 1, 2,

ALTERNATIVE - Comment on gait ALTHAROLD_042_001_010 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** false **Is Annotated:** false

Description: ALTERNATIVE - Comment_on_Gait

ALTERNATIVE - Pelvic elevation ALTHAROLD_042_001_011 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - Pelvic_elevation

Options: 0, 1, 2, 3,

ALTERNATIVE - Tail elevation ALTHAROLD_042_001_012 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - Tail_elevation

Options: 0, 1, 2,

ALTERNATIVE - Startle response ALTHAROLD_042_001_013 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - Startle_Response

Options: 0, 1, 2, 3,

ALTERNATIVE - Touch escape ALTHAROLD_042_001_014 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - Touch_Escape

Options: 0, 1, 2,

ALTERNATIVE - Positional passivity ALTHAROLD_042_001_015 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - Positional_Passivity

Options: 0, 1,

ALTERNATIVE - Trunk curl ALTHAROLD_042_001_016 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - Trunk_Curl

Options: 0, 1,

ALTERNATIVE - Limb grasping ALTHAROLD_042_001_017 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - Limb_Grasping

Options: 0, 1,

ALTERNATIVE - Pinna reflex ALTHAROLD_042_001_018 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** false **Is Annotated:** false

Description: ALTERNATIVE - Pinna_Reflex

Options: 0, 1,

ALTERNATIVE - Corneal reflex ALTHAROLD_042_001_019 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** false **Is Annotated:** false

Description: ALTERNATIVE - Corneal_Reflex

Options: 0, 1,

ALTERNATIVE - Contact righting reflex ALTHAROLD_042_001_020 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** false **Is Annotated:** false

Description: ALTERNATIVE - Contact_Righting_Reflex

Options: 0, 1,

ALTERNATIVE - Evidence of biting ALTHAROLD_042_001_021 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - Evidence_of_Biting

Options: 0, 1,

ALTERNATIVE - Vocalization ALTHAROLD_042_001_022 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - Vocalization

Options: 0, 1,

ALTERNATIVE - General comments about the mouse ALTHAROLD_042_001_023 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** false **Is Annotated:** false

Description: ALTERNATIVE - General_comments_about_the_mouse

ALTERNATIVE - Spontaneous activity ALTHAROLD_042_001_024 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - spontaneous_activity

Options: 0, 1, 2, 3,

ALTERNATIVE - Visual placing ALTHAROLD_042_001_025 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - visual_placing

Options: 0, 1, 2,

ALTERNATIVE - Grip strength ALTHAROLD_042_001_026 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - grip_strength

Options: 0, 1, 2, 3,

ALTERNATIVE - Body tone ALTHAROLD_042_001_027 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - body_tone

Options: 0, 1, 2,

ALTERNATIVE - Wire manoeuvre ALTHAROLD_042_001_028 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - wire_manoeuvre

Options: 0, 1, 2, 3, 4,

ALTERNATIVE - Heart rate ALTHAROLD_042_001_029 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - heart_rate

Options: 0, 1, 2,

ALTERNATIVE - Limb tone ALTHAROLD_042_001_030 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - limb_tone

Options: 0, 1, 2, 3,

ALTERNATIVE - Abdominal tone ALTHAROLD_042_001_031 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - abdominal_tone

Options: 0, 1, 2,

ALTERNATIVE - Rightling reflex ALTHAROLD_042_001_032 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - rightling_reflex

Options: 0, 1, 2, 3,

ALTERNATIVE - Catelepsy ALTHAROLD_042_001_033 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - catelepsy__geotaxis_

Options: 0, 1, 2,

ALTERNATIVE - Fear ALTHAROLD_042_001_034 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - fear

Options: 0, 1,

ALTERNATIVE - Irritability ALTHAROLD_042_001_035 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - irritability

Options: 0, 1,

ALTERNATIVE - Aggression ALTHAROLD_042_001_036 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - aggression

Options: 0, 1,

ALTERNATIVE - Date of procedure ALTHAROLD_042_001_037 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - date_of_procedure

ALTERNATIVE - Respiration rate ALTHAROLD_042_001_038 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - respiration_rate

Options: 0, 1, 2, 3,

ALTERNATIVE - Number of animals per cage ALTHAROLD_042_001_801 | v1.0

procedureMetadata

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Unit Measured: number

Description: ALTERNATIVE - Number_of_animals_per_cage

ALTERNATIVE - Location of test ALTHAROLD_042_001_802 | v1.0

procedureMetadata

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: ALTERNATIVE - Location_of_test

Options: laminar flow, open bench,

ALTERNATIVE - Number of days since animals were cleaned ALTHAROLD_042_001_803 | v1.0

procedureMetadata

Req. Analysis: false **Req. Upload:** false **Is Annotated:** false

Description: ALTERNATIVE - Number_of_days_since_animals_were_cleaned

---