Indirect ophthalmoscopy HAROLD_043_001

Parameters and Metadata

Eye HAROLD_043_001_001 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** true

Description: Eye

Options: normal, abnormal,

Description HAROLD_043_001_002 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** false **Is Annotated:** false

Description: Description

Left retina HAROLD_043_001_003 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** false **Is Annotated:** true

Description: Left_Retina

Options: normal, abnormal, no data,

Left retina pigmentation HAROLD_043_001_004 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** false **Is Annotated:** true

Description: Left_Retina_Pigmentation

Options: normal, abnormal, not defined, no data,

Left retina structure HAROLD_043_001_005 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** false **Is Annotated:** true

Description: Left_Retina_Structure

Options: normal, degenerate, not defined, no data,

Left optic disc HAROLD_043_001_006 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** false **Is Annotated:** true

Description: Left_Optic_Disc

Options: normal, abnormal, not defined, no data,

Right retina HAROLD_043_001_007 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** false **Is Annotated:** true

Description: Right_Retina

Options: normal, abnormal, no data,

Right retina pigmentation HAROLD_043_001_008 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** false **Is Annotated:** true

Description: Right_Retina_Pigmentation

Options: normal, abnormal, not defined, no data,

Right retina structure HAROLD_043_001_009 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** false **Is Annotated:** true

Description: Right_Retina_Structure

Options: normal, degenerate, not defined, no data,

Right optic disc HAROLD_043_001_010 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** false **Is Annotated:** true

Description: Right_Optic_Disc

Options: normal, abnormal, not defined, no data,

Left eye blood vessels HAROLD_043_001_011 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** false **Is Annotated:** true

Description: Left_Eye_Blood_Vessels

Options: normal, abnormal, no data,

Left eye blood vessels structure HAROLD_043_001_012 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** false **Is Annotated:** true

Description: Left_Eye_Blood_Vessels_Structure

Options: normal appearance, abnormal appearance, not defined, no data,

---

Left eye blood vessels pattern HAROLD_043_001_013 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** false **Is Annotated:** true

Description: Left_Eye_Blood_Vessels_Pattern

Options: normal, abnormal, not defined, no data,

Right eye blood vessels HAROLD_043_001_014 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** false **Is Annotated:** true

Description: Right_Eye_Blood_Vessels

Options: normal, abnormal, no data,

Right eye blood vessels structure HAROLD_043_001_015 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** false **Is Annotated:** true

Description: Right_Eye_Blood_Vessels_Structure

Options: normal appearance, abnormal appearance, not defined, no data,

---

Right eye blood vessels pattern HAROLD_043_001_016 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** false **Is Annotated:** true

Description: Right_Eye_Blood_Vessels_Pattern

Options: normal, abnormal, not defined, no data,

Digitial image(s) HAROLD_043_001_017 | v1.0

seriesMediaParameter

Req. Analysis: false **Req. Upload:** false **Is Annotated:** false

Description: Digitial_Image_s_

Body weight HAROLD_043_001_018 | v1.0

simpleParameter

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Unit Measured: g

Description: Body_Weight

Equipment name HAROLD_043_001_801 | v1.0

procedureMetadata

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: Equipment_name

Equipment manufacturer HAROLD_043_001_802 | v1.0

procedureMetadata

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: Equipment_manufacturer

Equipment model HAROLD_043_001_803 | v1.0

procedureMetadata

Req. Analysis: false **Req. Upload:** true **Is Annotated:** false

Description: Equipment_model

Topical agent(s) HAROLD_043_001_804 | v1.0

procedureMetadata

Req. Analysis: true **Req. Upload:** true **Is Annotated:** false

Description: Topical_agent_s_

Options: atropine, tropicamide, atropine + neosynephrine, tropicamide + neosynephrine, tropicamide + phenylephrine,

---